

# Gram staining technique

- 1) For smears from solid medium (e.g. nutrient agar plates or slopes, emulsify a very small portion of a single colony in a drop of 0.9% saline on a clean plain microscope slide. For smears from liquid media (e.g. nutrient broth), transfer a loopful of fluid to a clean microscope slide. This will contain far fewer organisms than smears of colonies from solid media.
- 2) Allow material on slide to completely air dry.
- 3) Pass the slide quickly through a Bunsen flame (approximately 0.5 seconds) to fix the material to the slide. Do not overheat.
- 4) Place a large drop of crystal violet on the slide, over the smear. Leave the stain on the slide for 1 minute.
- 5) Rinse the slide under running tap water for 5 seconds.

NOTE: At each rinsing, take care not to wash smear off the slide.

- 6) Cover the smear with Gram's iodine solution and leave for one minute.
- 7) Rinse the slide under running water for 5 seconds.
- 8) Decolorise briefly with acetone alcohol for 2 seconds.
- 9) Rinse under running water for 5 seconds.
- 10) Counter stain by completely covering the smear with Safranin stain for 30 seconds.
- 11) Rinse briefly with running water and blot slide dry. A piece of white clean paper can be used to blot the slide.
- 12) Gram positive species will appear blue to purple, while Gram negative species will appear red to pink.

Disposal: At the end of a microbiology activity it is important that all equipment be sterilised prior to disposal. Sterilise at 121°C and a pressure of 15 psi for 20-30 minutes prior to disposal. Verify conditions with a sterile confirmation strip (E5.36).

For safe laboratory handling and disposal, refer to AS/NZS2243 Safety in Laboratory Standards.