

## miniPCR Learning Lab™: PTC Taster Lab From Genotype to Phenotype

# The DNA Times Free

## Parents Genetically Engineer Children to Like Broccoli

**A diagnostic clinic in California is developing an assay to test couples for sensitivity to bitter chemicals as the first step to genetically engineering their children to like broccoli**

Boston, October 30<sup>th</sup>, 2037. Parents are tired of having to force feed their kids broccoli. Doctors are now offering the option of genetically engineering their unborn children to decrease their sensitivity to broccoli's bitter taste. After identifying the gene responsible for aversion to the bitter tasting compound found in broccoli, a new gene therapy is about to hit the market to correct the aversion in unborn children.



The Taqdale Clinic in California is developing a screening assay to test future parents to determine if they are candidates for gene therapy. The clinic is using phenylthiocarbamide (PTC), a bitter compound similar to that found in broccoli, that can only be tasted by people carrying specific variants of the *TAS2R38* gene. Student volunteers have been invited to help test the assay before its market release.

Woolly mammoth finally cloned at  
Brooklyn Zoo p.5



Transgenic heart-shaped bananas  
could solve world hunger p.7



## miniPCR Learning Lab™: PTC Taster Lab From Genotype to Phenotype

### Instructor's Guide Contents

1.	<a href="#"><u>Synopsis</u></a>	p. 3
2.	<a href="#"><u>Learning goals and skills developed</u></a>	p. 4
3.	<a href="#"><u>Standards alignment</u></a>	p. 5
4.	<a href="#"><u>Background and significance</u></a>	p. 7
5.	<a href="#"><u>Laboratory set-up manual</u></a>	p. 11
6.	<a href="#"><u>Laboratory guide</u></a>	p. 15
7.	<a href="#"><u>Study questions</u></a>	p. 26
8.	<a href="#"><u>Student-centered investigations and extension activities</u></a>	p. 31
9.	<a href="#"><u>Ordering information</u></a>	p. 32
10.	<a href="#"><u>About miniPCR Learning Labs™</u></a>	p. 33

*To taste or not to taste*, that is the question. And the answer is surprisingly simple: a few nucleotides. In this lab, we will examine how single nucleotide polymorphisms (SNPs) can change our ability to perceive the world around us. The human sense of taste is composed of an intricate neurophysiological network, but it only takes small changes to one gene to change the way we taste...

## 1. Synopsis

Approximately 99.9% of the human genome is identical across individuals; the other 0.1% makes each one of us unique. The human genome contains about 3 billion base pairs but just a single nucleotide change can alter or abolish gene function and can result in a new phenotype. A *single nucleotide polymorphism*, or SNP (pronounced "snip"), is a common variant at a single position in a DNA sequence among individuals of the same species. SNPs are the most common form of variation in the human genome and account for much of our genetic and phenotypic diversity.

In this Learning Lab, students will explore SNPs associated with their own phenotypes. Students will assess their ability to taste the chemical phenylthiocarbamide (PTC) and determine how that ability correlates with their genotype at the *TAS2R38* locus, which encodes for a taste receptor expressed in gustatory papillae. There are two common alleles for the *TAS2R38* gene, a 'taster' allele and a 'non-taster' allele. The difference between these alleles results from the combination of just three SNPs, or three single-nucleotide changes. It is a goal of this lab to illustrate how very small genetic changes can have significant functional consequences.

Students will extract their own DNA from cheek cells, perform PCR (polymerase chain reaction) to amplify the *TAS2R38* gene, use a restriction enzyme on the PCR product to distinguish between *TAS2R38* alleles, and perform gel electrophoresis to visualize the results. This activity allows students to directly explore the genetic basis of a fundamental aspect of neurophysiology—our ability to taste.

- **Techniques utilized:** DNA extraction, PCR, restriction endonuclease digestion, gel electrophoresis, and DNA visualization
- **Time required:** At least two 45-minute periods (see Section 6, *Planning your time*). This does not include PCR run time of approximately 90 minutes.
- **Reagents needed:** *miniPCR PTC Taster Lab* reagents kit (available at [miniPCR.com](http://miniPCR.com)), gel electrophoresis reagents (See sections 5 and 9)
- **Suggested skill level:** Familiarity with central dogma of molecular biology, molecular genetics, micropipetting techniques

## 2. Learning goals and skills developed

### Student Learning Goals – students will:

- Comprehend how DNA encodes traits that are transmitted across generations
- Understand that PCR is a technique for amplifying specific parts of the genome
- Learn about the molecular physiology of taste
- Learn about GPCRs signal transduction processes in neurophysiology
- Understand the link between genotype and phenotype
- Study the molecular basis of genetic inheritance
- Learn about the implications of SNPs in personalized medicine and in evolution

### Scientific Inquiry Skills – students should be able to:

- Formulate hypotheses and predict results
- Make predictions about their genotypes at a specific locus based on empirical observations of a trait (ability or inability to taste PTC)
- Compare genetic analysis results to their predictions and draw conclusions based on hypotheses
- Engage in an evidence-based debate around genetics and inheritance
- Generate data visualizations to present their results

### Molecular Biology Skills:

- Micropipetting
- Isolation of genomic DNA from cheek cells
- Polymerase Chain Reaction (PCR)
- Analysis of SNPs through restriction digest
- Preparation of agarose gels
- Agarose gel DNA electrophoresis
- Staining, visualization, and molecular weight analysis of DNA fragments

### 3. Standards alignment

**Next Generation Science Standards - Students will be able to...**

- HS-LS1-1 Construct an explanation based on evidence for how the structure of DNA determines the structure of proteins which carry out the essential functions of life through systems of specialized cells.
- HS-LS3-1 Ask questions to clarify relationships about the role of DNA and chromosomes in coding the instructions for characteristic traits passed from parents to offspring.
- HS-LS3-2 Make and defend a claim based on evidence about the natural world that reflects scientific knowledge, and student-generated evidence.
- HS-LS4-1 Communicate scientific information that common ancestry and biological evolution are supported by multiple lines of empirical evidence.
- HS-LS4-4 Construct an explanation based on evidence for how natural selection leads to adaptation of populations.

**Common Core English Language Arts Standards - Students will be able to...**

- RST.9-10.9 Compare and contrast findings presented in a text to those from other sources (including their own experiments), noting when the findings support or contradict previous explanations.
- RST.9-10.3 Follow precisely a complex multistep procedure when carrying out experiments, taking measurements, or performing technical tasks, attending to special cases or exceptions defined in the text.
- RST.11-12.1 Cite specific textual evidence to support analysis of science and technical texts, attending to important distinctions the author makes and to any gaps or inconsistencies in the account.
- WHST.9-12.2 Write informative/explanatory texts, including the narration of historical events, scientific procedures/ experiments, or technical processes.

WHST.9-12.9 Draw evidence from informational texts to support analysis, reflection, and research.

SL.11-12.5 Make strategic use of digital media (e.g., textual, graphical, audio, visual, and interactive elements) in presentations to enhance understanding of findings, reasoning, and evidence and to add interest.


## Advanced Placement® Biology Science Practices and Learning Objectives

- SP 1 The student can use representations and models to communicate scientific phenomena and solve scientific problems.
- SP 3 The student can engage in scientific questioning to extend thinking or to guide investigations within the context of the AP® course.
- SP 5 The student can perform data analysis and evaluation of evidence.
- SP 6 The student can work with scientific explanations and theories.
- SP 7 The student is able to connect and relate knowledge across various scales, concepts and representations in and across domains.
- LO 3.1 The student is able to construct scientific explanations that use the structures and mechanisms of DNA and RNA to support the claim that DNA and, in some cases, that RNA are the primary sources of heritable information.
- LO 3.6 The student can predict how a change in a specific DNA or RNA sequence can result in changes in gene expression.
- LO 3.16 The student is able to explain how the inheritance patterns of many traits cannot be accounted for by Mendelian genetics.
- LO 3.18 The student is able to describe the connection between the regulation of gene expression and observed differences between different kinds of organisms.
- LO 3.24 The student is able to predict how a change in genotype, when expressed as a phenotype, provides a variation that can be subject to natural selection.
- LO 3.25 The student can create a visual representation to illustrate how changes in a DNA nucleotide sequence can result in a change in the polypeptide produced.
- LO 3.26 The student is able to explain the connection between genetic variations in organisms and phenotypic variations in populations.
- LO 3.43 The student is able to construct an explanation, based on scientific theories and models, about how nervous systems detect external and internal signals, transmit and integrate information, and produce responses.

## 4. Background and significance

### The Molecular Physiology of Taste

G-protein-coupled receptors (GPCRs) are a large class of diverse and highly conserved proteins with seven transmembrane domains (Figure 1). Thousands of different human GPCRs have been identified. They play a role in a wide range of physiological processes, from hormonal communication to sensory perception. GPCRs are embedded in a cell's plasma membrane where they interact with G proteins that trigger signal transduction cascades inside the cell.



**Figure 1.** The seven-transmembrane structure of a G protein–coupled receptor, via [Wikipedia](#)

Our ability to taste is mediated by GPCRs. Taste receptor genes of the TAS1R and TAS2R families encode GPCRs expressed in taste receptor cells of the gustatory papillae in the tongue. Various taste receptors mediate specific taste modalities (bitter, sweet, savory, etc.). Humans have 43 different receptors for bitter taste within the TAS2R gene family. One of the best characterized taste receptors is encoded by the *TAS2R38* gene. This receptor binds ligands such as 6-n-propylthiouracil (PROP6) and phenylthiocarbamide (PTC), two chemicals similar to those found in bitter tasting vegetables. When ligands bind to the receptor they cause a conformational

change in the TAS2R38 receptor protein, which activates a signal transduction pathway that ultimately transmits a nerve impulse to the brain letting it know “that’s bitter”.

### To taste or not to taste... determined by a single gene

The TAS2R38 receptor is one of the only taste receptors for which variations at the DNA level are known to influence taste perception. Single nucleotide polymorphisms (SNPs) in the *TAS2R38* gene have been associated with impaired bitter taste perception (for a broad genomic analysis of bitter taste perception see Roudnitzky et al., 2015<sup>1</sup>). This is thought to be due to a deficit in signal transduction after the receptor binds to its ligands (such as the chemicals PROP6 and PTC). There are 3 SNPs commonly found on the *TAS2R38* gene that are correlated with an individual’s ability to taste PTC, resulting in two alleles often referred to as *PAV* and *AVI*, as detailed in Table 1.

Position (bp)	Position (amino acid)	Taster ( <i>PAV</i> )		Non-Taster ( <i>AVI</i> )	
		base	aa	base	aa
145	49	C	Pro	G	Ala
785	262	C	Ala	T	Val
886	296	G	Val	A	Ile

**Table 1.** Polymorphisms within the *TAS2R38* gene

The *PAV* allele is associated with high sensitivity to the bitter taste of PTC (‘Taster’ allele), while *AVI* is associated with little or no sensitivity (‘Non-taster’ allele). In human populations we find individuals that are homozygous for *PAV* (‘strong tasters’), homozygous for *AVI* (‘non-tasters’), or heterozygous (‘moderate tasters’).

### Implications of SNPs in human health

There are roughly 10 million SNPs in the human genome, and most of them have no known effect in health or development. Yet many SNPs have been linked to other forms of phenotypic

<sup>1</sup> Roudnitzky N, Behrens M, Engel A, Kohl S, Thalmann S, Hübner S, Lossow K, Wooding SP, Meyerhof W. Receptor Polymorphism and Genomic Structure Interact to Shape Bitter Taste Perception. PLoS Genet. 2015 Sep 25;11(9):e1005530. doi:10.1371/journal.pgen.1005530. eCollection 2015 Sep. PubMed PMID: 26406243; PubMed Central PMCID: PMC4583475.

variation besides sensitivity to PTC, including SNP association with certain diseases. In some cases SNPs may directly influence physiology (e.g. by directly altering the coding sequence for a gene), and in other cases they simply serve as biological markers that are associated with genetic traits (e.g. when they occur in the vicinity of a gene). Once phenotypic associations of SNPs are known, scientists can examine changes in the surrounding DNA in an attempt to identify the gene or genes responsible for the trait. Hence some SNPs can help predict a patient's risk of developing a certain disease.

Some SNPs can also determine a patient's differential response to a drug. About 50% of medicines on the market target GPCRs, therefore a better understanding of GPCRs and their associated polymorphisms may serve as a gateway to personalized medicine. Understanding a patient's genotype at a specific pharmacological target can help predict their response to treatment, enabling healthcare providers to harness our growing understanding of SNPs to inform their practice.

## Today's Lab

Students will determine the genotype of SNP785 of their own *TAS2R38* genes, in which a **cytosine to thymine** change at the nucleotide level results in an **alanine to valine** change at the amino acid level. At the same time students will test their ability to taste the bitter compound PTC and correlate it with their genotypes, which will later be tested through molecular analysis.

Students will obtain cheek cells by taking a swab of the inside of their mouths and will isolate their genomic DNA by lysing the cells. Using their genomic DNA as template, students will amplify a 250 bp DNA region surrounding SNP785 of *TAS2R38* using PCR. In order to distinguish between 'taster' and 'non-taster' alleles, students will perform a restriction digest of the PCR products using a restriction endonuclease that can only recognize and cleave the 'taster' allele (Figure 2). The cleaved DNA will result in two fragments of approximately 150 bp and 100 bp, while DNA corresponding to the 'non-taster' allele will remain intact at 250 bp. The difference in size will be visualized using gel electrophoresis.




Figure 2. PTC Taster Lab overview

## 5. Laboratory set-up manual

Reagent	Amount per group	Storage	Checklist
<b>miniPCR EZ PCR Master Mix, Load-Ready™</b> <ul style="list-style-type: none"> <li>• 2X Mix with <i>Taq</i> DNA polymerase</li> <li>• dNTPs (included)</li> <li>• PCR buffer with Mg<sup>2+</sup> (included)</li> <li>• Gel Loading Dye (included)</li> </ul>	60 µL	-20°C freezer	Supplied in Kit
<b>PTC Primer Mix</b> <ul style="list-style-type: none"> <li>• 2X concentrate</li> </ul>	60 µL	-20°C freezer	
<b>X-Tract™ Buffer</b> <ul style="list-style-type: none"> <li>• DNA extraction buffer, alkaline lysis</li> </ul>	240 µL	-20°C freezer or room temperature	
<b>100 bp DNA Ladder, Load-Ready™</b> <ul style="list-style-type: none"> <li>• Pre-mixed with gel loading dye</li> </ul>	15 µL	-20°C freezer	
<b>Restriction Enzyme</b> <ul style="list-style-type: none"> <li>• Fnu4HI</li> </ul>	5 µL	-20°C freezer	
<b>PTC Taste Paper</b> <ul style="list-style-type: none"> <li>• PTC strips</li> <li>• Control paper strips</li> </ul>	4 strips of each	Room temp.	
<b>Flat-Head Toothpicks</b>	4 sticks	Room temp.	
<b>DNA Staining Agent</b> <ul style="list-style-type: none"> <li>• e.g., GreenView Plus (for Blue light)</li> </ul>	Follow supplier instructions	4°C, dark	
<b>2% Agarose Gels</b> <ul style="list-style-type: none"> <li>• Electrophoresis grade agarose</li> </ul>	5 gel lanes per group (4 samples, plus 100 bp ladder)	Room temp.	
<b>1X TBE Electrophoresis Buffer</b>	50 ml per blueGel™	Room temp.	

Equipment and Supplies	Checklist
<b>PCR Thermal cycler:</b> e.g. miniPCR™ machine	
<b>Micropipettes</b> <ul style="list-style-type: none"> <li>• 0.5-10 µL: at least one for the classroom</li> <li>• 2-20 µL: one per lab group</li> <li>• 20-200 µL: one for the teacher</li> <li>• 100-1000 µL: one for the teacher</li> </ul>	
<b>Disposable micropipette tips</b>	
<b>PCR tubes:</b> 200 µL microtubes, individual or in 8-tube strips (12 tubes per group)	
<b>Plastic tubes:</b> 1.5 or 1.7 mL tubes to aliquot reagents (4 tubes per group)	
<b>Horizontal gel electrophoresis apparatus,</b> e.g. blueGel™ by miniPCR	
<b>DC power supply</b> for electrophoresis apparatus	
<b>Transilluminator:</b> Blue light recommended (or UV with protective gear) <ul style="list-style-type: none"> <li>• Not needed with blueGel™</li> </ul>	
<b>Scale</b> for weighing agarose	
<b>250ml flasks or beakers</b> to dissolve agarose gel	
<b>Microwave or hot plate</b> to prepare agarose gel	
<b>Microcentrifuge</b> (optional, to briefly spin down kit components before use)	
<b>Gel documentation system</b> (optional, or use cell phone camera instead)	
<b>Other supplies:</b> <ul style="list-style-type: none"> <li>• Disposable laboratory gloves</li> <li>• Permanent marker</li> <li>• Small beaker for material disposal</li> </ul>	

## Planning your time

This lab has 6 stages, and is designed to run in **two 45-min** class periods, or a **single 3-h block**

### Experimental Period 1 (~45 min)

*Preparation: Dispense reagents and prepare equipment*

- 20 min

- A DNA extraction**
  - 15 min
- B PCR set up**
  - 10 min
- C PCR programming & monitoring**
  - 5 min PCR programming
  - 15 min PCR monitoring, discussion




**Stopping point.** miniPCR will run for **approximately 1.5h**  
Once PCR is completed, you can store PCR product in fridge (up to 1 week) or freezer (longer term)

### Experimental Period 2 (~45 min)

*Preparation: Pour agarose gels*

- 20 min

- D Restriction digest, PTC Taste test**
  - 15 min
- E Gel electrophoresis**
  - 25 min
- F Size determination & interpretation**
  - 5 min



## Quick guide: Preparatory activities

Suggested for 8 Student Groups (4 Students per group), 32 total samples tested

### A. DNA Extraction

- Provide each group with 4 flat-end toothpicks or sterile loops
- Each group will need 4 PCR tubes (200  $\mu$ l)
- Aliquot 240  $\mu$ L of X-Tract™ DNA Extraction Buffer per group (or  $\mu$ L 50 per student)

### B. PCR set up

- Thaw primers and master mix by placing tubes on a rack or water bath at room temperature
- Label and dispense the following tubes per group (for 4 PCR reactions):
 

- X-Tract™ DNA Extraction Buffer	240 $\mu$ L
- 2X EZ PCR Master Mix	60 $\mu$ L
- 2X PTC Primer Mix	60 $\mu$ L
- Each lab group will additionally need:
  - One 2-20 $\mu$ L micropipette
  - Micropipette tips and a small beaker or cup to dispose them
  - 4 PCR tubes (thin walled, 200  $\mu$ L microtubes)
  - Permanent marker (fine-tipped)

### B. PCR programming and monitoring

- Ensure lab benches are set up with a miniPCR™ machine and power supply

### C. Restriction digest

- Each group will need 4 PCR tubes (thin walled, 200  $\mu$ L microtubes)
- One 0.5-10  $\mu$ L micropipette to dispense restriction enzyme (may be shared by the class)

### E. PTC taste test

- Each student will need 1 PTC paper strip
- Print out a table or spreadsheet to record taster phenotypes (e.g., non-taster) and genotypes

### F. Gel electrophoresis

- Gels can be poured in advance of the class (described in section 6E)
- 5 Lanes per group (4 Test Samples + 1 lane for DNA ladder), e.g. 40 lanes for 32 students
- Pre-poured gels can be stored in the fridge, in a sealed container or in plastic wrap, in the dark
- If doing the gel run on a different day than the PCR run, completed PCR tubes can be stored in the fridge for up to one week until they are used, or in the freezer for longer-term storage

### G. Size determination and interpretation

- Have the banding pattern of the 100bp DNA Ladder handy to help interpret the electrophoresis results (see page 23)

## 6. Laboratory guide

### A. DNA Extraction

#### 1. Extraction Tube Preparation:

- Dispense **1 PCR tube (200  $\mu$ L thin-walled PCR tube) per student**
- **Label tube on the side, not cap, of the tube** with the student's initials

#### 2. Add 50 $\mu$ L of X-Tract™ DNA Extraction Buffer to each 200 $\mu$ L tube



**NOTE: for maximum DNA extraction efficiency, please abstain from eating or chewing gum ~20 minutes prior to cheek cell collection**

#### 3. Scrape the inside of your cheek 5-6 times with a flat-end toothpick to collect cells at the end of the toothpick



Rub *gently* along cheek, it shouldn't hurt!

#### 4. Dip the toothpick in the 200 $\mu$ L tube swirling it in X-Tract™ Buffer

- Swirl toothpick thoroughly in the buffer to release cells

#### 5. Tightly cap the tubes

#### 6. Incubate the tube in miniPCR™ for 10 minutes at 95°C

- Insert tubes and close the PCR machine lid and gently tighten the lid using the adjustment knob
- Use miniPCR™ Heat Block mode, or alternatively use a 95°C heat block or water bath. *If you use a heat block or water bath make sure to put pressure on the lids, otherwise the lids might pop open.*

#### 7. Remove tubes from heat and let them rest in a tube rack at room temperature



We recommend using the DNA extract immediately for PCR

## B. PCR set up

1. **PCR Tube Preparation:** Label 4 clean PCR tubes (200  $\mu$ L thin-walled tubes) per group
  - Student's initials on the side wall of the tube



*Remember to change tips at each of the following steps!*

2. **Add PCR reagents** to each tube
  - 12.5  $\mu$ L PTC Primer Mix
  - 12.5  $\mu$ L EZ PCR Master Mix

PCR Tube (Test Groups)	
PTC Primer Mix	12.5 $\mu$ L
2X EZ PCR Master Mix	12.5 $\mu$ L
Template DNA	3 $\mu$ L
	28 $\mu$ L

3. **Add DNA samples** to each new tube

- 3 $\mu$ L Student DNA extracts from the top of the solution
- Pipette up and down to mix directly into the PCR solution



If there is a cloudy precipitate at the bottom of the tube *be careful not to transfer it to the PCR set up*

4. **Cap the tubes**


- Make sure all the liquid volume collects at the bottom of the tube
- If necessary, spin the tubes briefly using a microcentrifuge, or tap them on the benchtop
- Press firmly on the tube caps to ensure a tight fit

5. **Place the tubes inside the miniPCR™ machine**

- Close the PCR machine lid and gently tighten the lid

## PCR programming and monitoring (illustrated using miniPCR™ software)


1. Open the miniPCR software app and remain on the "Protocol Library" tab



2. Click the "New Protocol" button on the lower left corner
3. Select the PCR "Protocol Type" from the top drop-down menu
4. Enter a name for the Protocol; for example "PTC Taster Lab – Group 1"
5. Enter the PCR protocol parameters:

- **Initial Denaturation**                    94°C    120 sec
- **Denaturation**                            94°C    10 sec
- **Annealing**                                58°C    15 sec
- **Extension**                                 72°C    40 sec
- **Number of Cycles**                       30
- **Final Extension**                         72°C    5 sec
- **Heated Lid**                                ON

6. Click **"Save"** to store the protocol
7. Click **"Upload to miniPCR"** (and select the name of your miniPCR machine in the pop-up dialogue window) to finish programming the thermal cycler.
8. Make sure that the **power switch in the back** of miniPCR is in the **ON** position
9. Click on **"miniPCR [machine name]"** tab to begin monitoring the PCR reaction (the miniPCR serial number can be found on the back of the miniPCR).



*The miniPCR™ software allows each lab group to monitor the reaction parameters in real time, and to export the reaction data for analysis as a spreadsheet.*

10. Once the PCR run is completed (approximately 90 min), the screen will display: **"Status: Completed"**. All LEDs on the miniPCR machine will light up.



You can now open the miniPCR lid and remove your tubes.  
**Be careful not to touch the metal lid which may still be hot**

## C. Restriction Digest and PTC Taste Test

1. **Digest Tube Preparation: Label 4 clean PCR tubes** (200  $\mu$ L thin-walled tubes) per group
  - Student's initials on the side wall, followed by "-R" (e.g. SK-R)
2. **Add 14  $\mu$ L of PCR Product** to the bottom of each Digest Tube
3. **Add 1  $\mu$ L of Restriction Enzyme** to each Digest Tube, directly into the PCR product.
  - Use a 0.5-10.0  $\mu$ L range micropipette



Pipette up and down to mix well. This mixing is critical for efficiency of the digest

	Volume
PCR Product	14 $\mu$ L
Restriction Enzyme	1 $\mu$ L
<b>FINAL VOLUME</b>	<b>15 <math>\mu</math>L</b>

4. **Cap the tubes**
  - Make sure all the liquid volume collects at the bottom of the tube
  - If necessary, spin the tubes briefly using a microcentrifuge
5. **Place the tubes inside the miniPCR machine**
  - Press firmly on the tube caps to ensure a tight fit
  - Close the PCR machine lid and gently tighten the lid
6. Incubate in the miniPCR machine **15 minutes at 37°C**
  - a. Open the miniPCR software app and remain on the "**Protocol Library**" tab
  - b. Click the "**New Protocol**" button on the lower left corner
  - c. Select the "**Heat Block**" **Protocol Type** from the top drop-down menu
  - d. Enter a name for the Protocol; for example "**Restriction Digest – Group 1**"
  - e. Enter the Heat Block parameters:
    - **Heat Block**                      37°C, 15 minutes
  - f. Click "**Save**" to store the protocol
  - g. Click "**Upload to miniPCR**" (and select the name of your miniPCR machine in the dialogue window) to finish programming
  - h. Make sure that the **power switch in the back** of miniPCR is in the **ON** position
  - i. Click on "**miniPCR [machine name]**" tab to begin monitoring the incubation

## 7. ***PTC taste test (while Restriction digest is running)***

- a. Begin by asking students to taste “control” paper strips to establish a baseline, this can be any chemical-free paper (e.g. thin strips of standard printer paper)
- b. Pass out PTC paper strips to students, one per student
- c. Ask students to taste paper, and record their ability to taste the PTC chemical (i.e. *non taster, weak taster, strong taster*)
- d. Collect class data, if possible

Student Name	Assessed Phenotype	Predict Genotype	Actual Genotype
	Strong taster		
	Moderate taster		
	Non-taster		
	.....		

8. Once the Restriction Digest is completed (15 min), open miniPCR lid and remove tubes
9. The restriction digest must be used immediately for Gel Electrophoresis



Because Fnu4HI enzyme cannot be heat inactivated, *we do not recommend storing the digest for use at a later date*, as it may result in sample degradation

## D. Gel electrophoresis – Pouring agarose gels (Preparatory activity)

*If this lab is going to be performed over two class periods, gels can be prepared up to one day ahead of the second period and stored in a refrigerator, covered in plastic wrap and protected from light.*

1. Prepare a clean and dry agarose gel casting tray
  - Seal off the ends of the tray as indicated for your apparatus (not needed for blueGel™ users)
  - Place a well-forming comb at the top of the gel (5 lanes or more)
2. For each lab group, prepare a 2% agarose gel using electrophoresis buffer
  - Adjust volumes and weights according to the size of your gel tray
  - For example, add 1 g of agarose to 50 ml of electrophoresis buffer
  - **blueGel™ users:** add 0.4 gr of agarose to 20 ml of 1X TBE buffer
  - Mix reagents in glass flask or beaker and swirl to mix
  - If pouring more than one gel you can dissolve all the agarose at once (e.g. 4g in 200 ml for 10 blueGels)
3. Heat the mixture using a microwave or hot plate
  - Until agarose powder is dissolved and the solution becomes clear
  - Use caution, as the mix tends to bubble over the top and is very hot
4. Let the agarose solution cool for about 2-3 min at room temperature.
  - Swirl the flask intermittently
5. Add gel staining dye (e.g. GreenView™ Plus available at miniPCR.com)
  - Follow dye manufacturer instructions
  - Typically, 10 µL of staining dye per 100 mL of agarose solution
  - Note: We recommend the use of safe alternatives to ethidium bromide.
6. Pour the agarose solution into the gel-casting tray with comb
7. Allow gel to completely solidify (until firm to the touch) and remove the comb
  - Typically, ~10 minutes
8. Place the gel into the electrophoresis chamber and cover it with 1X TBE electrophoresis buffer

## ***Gel electrophoresis – Running the gel***

1. Make sure the gel is completely submerged in electrophoresis buffer
    - Ensure that there are no air bubbles in the wells (shake the gel gently if bubbles need to be dislodged)
    - Fill all reservoirs of the electrophoresis chamber and add just enough buffer to cover the gel and wells
    - **blueGel™ users:** Use no more than 25 ml of 1X TBE electrophoresis buffer
  
  2. Load DNA samples onto the gel in the following sequence
    - **Lane 1:** 10µL DNA ladder
    - **Lane 2:** 12µL Restriction Digest from **Student Sample 1**
    - **Lane 3:** 12µL Restriction Digest from **Student Sample 2**
    - **Lane 4:** 12µL Restriction Digest from **Student Sample 3**
    - **Lane 5:** 12µL Restriction Digest from **Student Sample 4**
- \* You will need 5 Lanes per group (i.e. 1 Ladder + 4 Samples)


Note: there is **no need to add gel loading dye to your samples.**

The *miniPCR EZ PCR Master Mix* and *100 bp DNA Ladder Load-Ready™* come premixed with loading dye, and ready to load onto your gel!

3. Place the cover on the gel electrophoresis box
  - Ensure the electrode terminals fit snugly into place
  
4. Set the voltage at 100-130V and conduct electrophoresis for 15-20 minutes, or until the colored dye has progressed to at about half the length of the gel (if using blueGel™, simply press the “Run” button.)
  - Check that small bubbles are forming near the terminals in the box
  - Longer electrophoresis times will result in better size resolution
  
5. Once electrophoresis is completed, turn the power off and remove the gel from the box

## E. Size determination and interpretation

1. Place the gel on the blue-light transilluminator
  - If using blueGel™ simply press the illuminator button
  - *If using UV light, cover any exposed skin and cover eyes with UV-protective goggles*
2. Verify the presence of PCR product
3. Ensure there is sufficient DNA band resolution in the 100-300 base pair (bp) range of the 100bp DNA ladder
  - Run the gel longer if needed to increase resolution
  - DNA ladder should look approximately as shown



*100 bp DNA Ladder visualized by ethidium bromide staining on a 1.3% TAE agarose gel.*

*Source: New England Biolabs*

4. Document the size of the PCR amplified DNA fragments by comparing the PCR products to the molecular weight reference marker (100 bp DNA ladder)
  - Capture an image with a smartphone camera or if available, use a Gel Documentation system

### F. Expected results and troubleshooting guide


Before restriction digest, PCR product for both the *taster* and *non-taster* alleles should migrate on agarose gel electrophoresis as a **single band at 250 base pairs (bp)**.

After PCR product has been cut with the restriction enzyme Fnu4HI:

- **The *non-taster allele* should migrate at 250bp**
- The *taster allele* should migrate as two separate fragments,
  - **one 150bp band, and one 100bp band**

Hence, there will be two common outcomes for *phenotypic tasters*:

1. A *homozygous taster* will typically display **two bands**: one at 150bp and one at 100bp
2. A *heterozygous taster* will display **three bands**: one at 250bp, one at 150bp, and one at 100bp
3. A *non-taster* will display only **one band**: 250bp



Presence of spurious, low-molecular weight bands at the bottom of the gel is not uncommon. These bands are usually “primer dimers” resulting from PCR primers binding onto each other and not the template DNA, replicating the primer itself and not the target DNA sequence.

### Troubleshooting

**The absence of a band may be attributed to inefficiency at different steps:**

- The **DNA was not extracted thoroughly**—no cells on the toothpick
- The **cheek cells were not properly dislodged** into the X-tract Buffer
- **Inefficient alkaline cell lysis** and DNA extraction
- **Insufficient PCR amplification**: suboptimal primer, Taq polymerase, nucleotides, or PCR parameters
- **The restriction digest was incomplete**: not long enough, insufficient enzymatic activity of the restriction endonuclease (Fnu4HI)
  - A faint band at 250 after digestion along with the presence of the 150 and 100bp bands **COULD** be the result of an incomplete digest, **NOT** a heterozygous genotype.
  - On the other hand, excessive incubation with some restriction enzymes may lead to sample degradation due to non-specific nuclease activity (star activity)
- Insufficient electrophoresis run or gel resolution
  - A thick band (double band) at ~150 bp after the digest could represent both the 150 bp and 100 bp bands: Let the gel run longer for better separation, or run again at higher agarose concentration (2.4%)
- Other sources of experimental error
  - Describe them:

## 7. Study questions - Pre-lab

### 1. Genetics of PTC sensitivity

PTC, or **1.** \_\_\_\_\_, was discovered accidentally in a lab by the scientists Arthur L. Fox and C.R. Noller. One night in his lab Fox accidentally over-poured PTC and complained of the **2.** \_\_\_\_\_ taste in the air as the powder diffused around him. His colleague Noller, on the other hand, tasted **3.** \_\_\_\_\_. Intrigued by the vast differences in their tasting Fox and Noller went on to explore the phenomenon further and to find an **4.** \_\_\_\_\_ for why only one of them could taste PTC while the other could not. They then gathered hundreds of people and had them sample PTC. They found that about 25% of respondents said they did not taste anything. Through further experimentation they ultimately found that a single gene was responsible for the **5.** \_\_\_\_\_ in taste perception of PTC. We now know that polymorphisms in the *TAS2R38* receptor are associated with differences in bitter taste perception. This phenotypic difference is passed down from **6.** \_\_\_\_\_ to **7.** \_\_\_\_\_ in a Mendelian fashion. At the DNA level, the two most common **8.** \_\_\_\_\_ for PTC differ in three **9.** \_\_\_\_\_. Each of these nucleotide differences is known as a SNP, for **10.** \_\_\_\_\_ .

### 2. Mendelian inheritance of PTC sensitivity

Sensitivity to the bitter taste of PTC is passed along from parents to offspring in a Mendelian fashion. This pattern can be analyzed using a Punnett square. In the example below, we observe the four possible genotypes for the offspring of two parents for the "Taster Gene", where "T" represents the dominant allele (Taster phenotype) and "t" represents the recessive allele (Non-Taster).

Suppose a heterozygous parent marries someone who is also heterozygous at the "Taster" locus.

- What genotypes would you expect in the progeny? What phenotypes?
- What type of cross is this? What allelic frequencies would you expect?

Punnett Square		Parent 1	
		Gametes	T
Parent 2	T		
	t		

### 3. Neurophysiology of taste

Taste receptors, or gustatory cells, are found in the **1.**\_\_\_\_\_ clustered in a group called a **2.**\_\_\_\_\_ \_\_\_\_\_. Gustatory papillae convey to the brain whether a food is salty or sweet, among other qualities. Humans have about **3.**\_\_\_\_\_ different receptors for bitter taste. These receptors are transmembrane proteins of the **4.**\_\_\_\_-\_\_\_\_\_ coupled receptor family.

### 4. Predict your taster phenotype and genotype

- a. Predict whether or not you will exhibit the PTC taster phenotype.
- b. If you are a taster of PTC, what are your possible genotypes at the *TAS2R38* locus?
- c. In which ways can single nucleotide polymorphisms (SNPs) affect the function of a gene?
- d. What do we mean when we say a SNP is *associated* with a certain phenotypic trait?
- e. Why do you think the three SNPs in the *TAS2R38* gene are inherited together (haplotype)?

### 5. G-Protein Coupled Receptors (GPCRs)

- a. What is the subcellular localization of GPCRs?
- b. List 3 known functions of GPCRs.
- c. List 3 types of cells that might be rich in GPCRs.

## Study Questions - During the Lab

### 1. DNA Extraction

- a. During DNA extraction we heated the cheek cells to 95°C. Why?
- b. Explain the role of the following components of a DNA extraction buffer:
  - **Detergent:**
  - **Chelating agent:**
  - **Basic (alkaline) ions:**

### 2. Polymerase Chain Reaction (PCR)

- c. What is happening to DNA molecules at each step of the PCR cycle?
  - **Denaturation:**
  - **Annealing:**
  - **Extension:**
- d. What is the purpose of the following PCR mix components?
  - **dNTPs:**
  - **Mg<sup>+2</sup>:**
  - **Taq DNA polymerase:**
- e. What temperature do you think is optimal for most enzymes?

- What feature of Taq DNA polymerase makes this enzyme uniquely suited for PCR?
  
- f. How many additional DNA molecules are there at the end of each PCR cycle? If the PCR consists of 30 cycles, how much more DNA will there be after the PCR completes?
  
- g. Which of the following is NOT a characteristic of PCR primers?
  - A. Short synthetic oligonucleotide
  
  - B. Typically 18-25 bases in length
  
  - C. Double stranded DNA
  
  - D. Unique homology to the DNA template
  
  - E. Sequence with ~50% G:C content

### 3. Restriction Digest

- a. Shown below is a partial sequence of your expected PCR product for the taster *TAS2R38* allele. Draw a box around the Fnu4HI restriction enzyme recognition site. Draw lines to illustrate where you expect the restriction enzyme will cut. How would the non-taster allele differ?

```

5' C T G C T T C T T T G T G A T A T C A T C C T G T G C T G C C T T C A T C T C T G T G C C C C T A C T 3'
3' G A C G A A G A A A C A C T A T A G T A G G A C A C G A C G G A A G T A G A G A C A C G G G G A T G A 5'
  
```

- b. What type of chemical bond is the Restriction Enzyme hydrolyzing (breaking)?
  
- c. Why is this enzyme called a *restriction endonuclease*? How is it different from an *exonuclease*?
  
- d. Why do we keep enzymes on ice when not in use? What does ice do to enzymatic activity?
  
- e. After the restriction digest, what is the expected size of the DNA fragments for the "Taster" allele?

## Study questions - Post Lab

- a. After tasting PTC paper, what did you expect your Taster genotype to be?
- b. Did this match your experimental results?
- c. What percentage of your class could taste PTC? What percentage couldn't? Did the experimental results match these distributions?
- d. What would be the expected allelic frequencies in a T/t x T/t cross? (Taster/non-Taster)
- e. Speculate what adaptive advantage might exist for animals able to taste PTC
- f. Speculate what adaptive implications might exist for humans able or unable to taste PTC
- g. Do you think the allelic frequencies of Taster and non-Taster vary across populations? Why?
- h. Why might someone with a heterozygous genotype find that they are unable to taste PTC?

## 8. Student-centered investigations and extension activities

### Student Centered Investigations

**Investigation 1:** Calculate Hardy-Weinberg equilibrium for your class data.

- Aggregating class data will give you counts of each genotype.
- Use genotype data to count the total number of each allele in your student population. This can be used to calculate the proportion (p and q) for each.
- Use the allele data to calculate the predicted genotype distribution using the Hardy-Weinberg equation.
- Compare the predicted genotype distribution with the actual distribution using a  $\chi^2$  test.
- An excel spreadsheet is available [here](#) on the miniPCR website to use as a data table.

## 9. Ordering information

### To order miniPCR™ PTC Taster Lab Kits:

- Phone: (781)-990-8PCR
- Email: [orders@minipcr.com](mailto:orders@minipcr.com)
- Online: [www.minipcr.com](http://www.minipcr.com)

### miniPCR™ PTC Taster Lab Kit (catalog no. KT-1004-03) contains:

*Materials are sufficient for a class of 32 students working in 8 lab groups*

- 2X EZ PCR Master Mix, Load-Ready™
  - includes *Taq* DNA polymerase, dNTPs, PCR buffer, and gel-loading dye
- PTC Lab 2X Primer Mix
- Restriction Enzyme Fnu4HI
- 100bp DNA Ladder, Load-Ready™
- X-Tract™ Buffer (*can be stored at ambient temperature*)

*These components should be kept frozen at -20°C*

*Reagents must be used within 12 months of shipment*

### Other reagents needed (not included in the kit):

- Agarose (electrophoresis grade)
- DNA staining agent (e.g., GreenView Plus available at [minipcr.com](http://minipcr.com))
- Gel electrophoresis buffer (e.g., TBE, available as 20X stock at [minipcr.com](http://minipcr.com))
- Distilled or deionized H<sub>2</sub>O
- Toothpicks (flat-headed)
- PTC paper (available on [Amazon.com](http://Amazon.com))

## 10. About miniPCR™ Learning Labs

This Learning Lab was developed by miniPCR™ in an effort to make molecular biology and molecular physiology more accessible.

We believe that there is no replacement for hands-on experimentation in the science learning process. Our goal is to augment everyone's love of scientific inquiry and STEM.

We develop miniPCR Learning Labs™ to help achieve this goal, working closely with educators, students, academic researchers, and others committed to science education. The guiding premise for this lab protocol is that a PCR-based lab plan based on real-life human genetics provides the right balance between intellectual engagement, experimentation, inquiry, and discussion. The design of this lab has simplified certain elements to achieve these goals: the analysis of bitter-tasting TAS2R38 haplotypes is reduced to variability at the SNP785 locus<sup>2</sup>. We use restriction digest analysis of a single SNP rather than a more complex restriction analysis to simplify interpretation and discussion of the results. Similar designs have been incredibly effective for other very successful educational lab courses, to which we owe inspiration<sup>3</sup>.

Starting on a modest scale working with Massachusetts AP Biology Summer Institute educators, this miniPCR™ Learning Lab has been received well, and its use is growing rapidly through academic and outreach collaborations.

*Authors:* Ezequiel Alvarez Saavedra, Ph.D., Sebastian Kraves, Ph.D., Leslie Prudhomme, Mandi Nyambi.

---

<sup>2</sup> For a detailed discussion of the genetics of bitter taste perception, see: Roudnitzky N., *et al.*, Receptor Polymorphism and Genomic Structure Interact to Shape Bitter Taste Perception. *PLoS Genet.* 2015 Sep 25;11(9):e1005530. doi:10.1371/journal.pgen.1005530. eCollection 2015 Sep. PubMed PMID: 26406243; PubMed Central PMCID: PMC4583475.

<sup>3</sup> See, for example: Bouakaze C, *et al.* "OpenLAB": A 2-hour PCR-based practical for high school students. *Biochem Mol Biol Educ.* 2010 Sep;38(5):296-302. doi: 10.1002/bmb.20408. PubMed PMID: 21567848.